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- [1] YangLiu,Xinyu Ye,Hong Zhang, et al. Antimicrobial Susceptibility of Mycoplasma pneumoniae Isolates and Molecular Analysis of Macrolide-Resistant Strains fromShanghai, China. ANTIMICROBIAL AGENTS AND CHEMOTHERAPY.2009,53(5):2160–2162 .
- [2] Rours GI, Hammerschlag MR, Van Doornum GJ, et al.Chlamydia trachomatis respiratory infection in Dutch infants. Arch Dis Child. 2009,94(9):705-7.
- [3] Nasution TA, Cheong SF, Lim CT, Multiplex PCR for the detection of urogenital pathogens in mothers and newborns. Malays J Pathol. 2007,29(1):19-24.
- [4] DeliXin, Zuhuang Mi,Xu Han, et al. Molecular Mechanisms of Macrolide Resistance in Clinical Isolates of Mycoplasma pneumoniae from China. ANTIMICROBIAL AGENTS AND CHEMOTHERAPY.2009,53(5):2158–2159.
- [5] Khan ER, Hossain MA, Paul SK, et al. Molecular diagnosis of genital Chlamydia trachomatis infection by polymerase chain reaction. Mymensingh Med J. 2011,20(3):362-5.
- [6] Robin Brittain-Long, Sandra Nord, Sigvard Olofsson, et al. Multiplex real-time PCR for detection of respiratory tract infections. Journal of Clinical Virology,2008, 1:53–56.
- [7] Maass V, Kern JM, Poeckl M, et al.Sequence Homologies between Mycoplasma and Chlamydia spp. Lead to False-Positive Results in Chlamydial Cell Cultures Tested for Mycoplasma Contamination with a Commercial PCR Assay. J Clin Microbiol. 2011,49(10):3681-2.
- [8] Brittain-Long R, Andersson LM, Olofsson S, et al. Seasonal variations of 15 respiratory agents illustrated by the application of a multiplex polymerase chain reaction assay. Scand J Infect Dis. 2011 Aug 26. [Epub ahead of print]
- [9] Thurman KA, Warner AK, Cowart KC, et al. Detection of Mycoplasma pneumoniae, Chlamydia pneumoniae, and Legionella spp. in clinical specimens using a single-tube multiplex real-time PCR assay. Diagn Microbiol Infect Dis. 2011 ,70(1):1-9.
- [10] Mitchell CM, Hovis KM, Bavoil PM, et al.Comparison of koala LPCoLN and human strains of Chlamydia pneumoniae highlights extended genetic diversity in the species. BMC Genomics. 2010,21;11:442.

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