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- [1] Tan EL, Chow VT, Kumarasinghe G, et al. Specific detection of enterovirus 71 directly from clinical specimens using real-time RT-PCR hybridization probe assay. J Mol Cell Probes. 2006;20 (2):135-140.
- [2] Tan EL, Yong LL, Quak SH, et al. Rapid detection of enterovirus 71 by real-time TaqMan RT-PCR. J Clin Virol. 2008;42 (2):203-206.
- [3] Fujimoto T, Yoshida S, Munemura T, Taniguchi K, et al. Detection and quantification of enterovirus 71 genome from cerebrospinal fluid of an encephalitis patient by PCR applications. J Infect Dis. 2008;61 (6):497-499.
- [4] Xiao XL, He YQ, Yu YG, et al. Simultaneous detection of human enterovirus 71 and coxsackievirus A16 in clinical specimens by multiplex real-time PCR with an internal amplification control. J Arch Virol. 2009; 154 (1):121-125.
- [5] Leś K, Przybylski M, Dzieciatkowski T, et al. Detection of human enteroviruses with real-time PCR assay using TaqMan fluorescent probe. J Med Dosw Mikrobiol. 2010; 62(3):245-253.
- [6] Marian AC Piqueur, Walter A Verstrepen, Peggy Bruynseels, et al. Improvement of a real-time RT-PCR assay for the detection of enterovirus RNA. Virol J. 2009; 6: 95.
- [7] Corless CE, Guiver M, Borrow R, et al. Development and evaluation of a “real-time” RT-PCR for the detection of enterovirus and parechovirus RNA in CSF and throat swab samples. J Med Virol. 2002;67(4):555-562.
- [8] Dierksen U, Rehren F, Henke-Gendo C, et al. Rapid routine detection of enterovirus RNA in cerebrospinal fluid by a one-step real-time RT-PCR assay. J Clin Virol. 2008;42(1):58-64.
- [9] Verstrepen WA, Bruynseels P, Mertens AH. Evaluation of a rapid real-time RT-PCR assay for detection of enterovirus RNA in cerebrospinal fluid specimens. J Clin Virol. 2002;25

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